



Pressing mode Kneading mode

Main unit scroll wheel

On/off: press and hold for 2 seconds Intensity: scroll from level 0 to 9 Mode: Short press, 4 modes are available



Hot compress function key



On/off/level switch: short press Audio prompt on/off: press and hold for 2 seconds Three levels of temperature are available



How to wear it correctly Wear the scroll wheel side to the left of your neck

03 EN

2. Press and hold scroll wheel for 2 seconds and turn on the main unit

with voice prompt. The default level is 0. It increases by 1 level when

3. Short press the scroll wheel of the main unit to cycle massage modes:

pressing mode, kneading mode, beating mode, integrated mode.

4. Press and hold the scroll wheel for 2 seconds to turn off the device.

you scroll the wheel up; When you scroll the wheel down, it decreases

by 1 level. The levels vary from 0-9. During scrolling, the main unit will

For each switch of mode, the intensity returns to default level 0. When

it's turned on again, the device will restore the massage mode chosen

5. Short press the Hot compress key to cycle the heating levels: low heat,

The temperatures corresponding to low-level heating, medium-level

6. The main unit is timed for 15 minutes and will stop massage automatically

7. When the main unit is on, it will turn off automatically after 120 seconds

1. When the main unit's battery runs low, the orange light of power button

charging port of the main unit and plug the other side into mobile phone

charger (output 5V 1-2A). When charging, the orange light is normally on.

is normally on. Plug one side of MicroUSB charging cable into the

When it's fully charged, the blue light is normally on.

if the product is not properly worn or the main unit is not operated.

heating and high-level heating are 38 °C, 40 °C, 42 °C respectively.

medium heat, high heat, heating off. The main unit will give audio prompt

1. First, wet the skin of neck part and wear the product.

give voice prompt of current level.

last time before it's turned off.

after 15 minutes.

Product charging

of current heating level after you operate.

Main unit operation

- Appropriate massage is recommended; it is recommended to use no
 - Store the device in a well-ventilated and dry place after use.
 - Don't put other heavy items on the device.

- When cleaning and maintaining the device, turn it off and clean the massage pads with wet tissue paper. When it's especially dirty, use a small amount of medical alcohol (alcohol concentration of 75%) to wide
- During daily maintenance, use neutral cleanser and soft cloth to wipe. Don't use corrosive cleanser and don't soak or rinse it in water. Don't
- wash it with water. During cleaning process, prevent water or other liquid from entering the device to avoid short circuit and electric shock.

- If you have problems using the product, please refer to the following tips • If the product is still in a fault condition or cannot be used, turn off the
- Users should not disassemble or repair the product by themselves. Dismantling or repairing by themselves will cause damage to the product and cause the invalidity of warranty. The company will not be responsible for any consequences caused thereby.

 Whether the internal rechargeable battery runs low. Unable to start • Please confirm the product power switch is on. Failure to charge | Please check if the USB cable and host and phone or no charge light power adapter are connected. Does the pad contact the skin directly? Adjust the position so that it is close to the skin. Normal display | • Is there hair or clothing, jewelry, etc. between the pads and the skin? Move away or remove items. stimulation • Is the intensity level at 0? Increase the massage intensity. • Is the massage pads and neck skin closely contacted? Adjust the position so that it's stick to the skin. Is the massage time too long? It is recommended to massage for 15 minutes at a time, twice a day. Is the massage intensity too high? Adjust to a low level • Is the skin allergic to the massage pads? Those with severe allergies need to be desensitized first. • Is the climate cold and the skin dry? Use a wet towel to moisten the skin area, prior to using. • Does the skin have a wound? Do not plate the massage pads on the wound during use. • Is the massage pads in good contact with the skin area? Try adjusting the position of the massage pads onto the skin surface, to improve skin contact, and adjust the massager's intensity. Weak simulation • Is the massage pads dirty or oily? Clean the massage pads surface. • Is the battery running low? Please charge if the battery is low. Is massage intensity level too low? Increase massage intensity.

Possible cause/treatment

05 EN

This product is a household product. Not intended for commercial use.

As we continually improve our products, we shall reserve the right to make

Thank you for choosing our Smart Neck Massager. Please read this manual

carefully before use in order to use the product correctly and take full

and comfort to your life.

Product diagram...

Precautions.....

Warranty card

advantage of its excellent performance. We believe it will bring convenience

Please keep this manual for your future reference. If any problems occur

when using the product, please contact SKG after-sales service outlet or our

company. We will be glad to help you. Thank you!

Operation instruction.....

Electrical schematic diagram.....

technical changes without prior notice.

Product and component recycling instruction....

Routine maintenance......

Technical parameters......

Note for discarding product.....

Danger Please never use the product together with the following medical

* It may cause malfunction of these medical electronic devices,

electronic devices: 1. Pacemaker and other implantable medical electronic devices; 2. Medical electronic devices such as ECMO that sustain life; 3. Electrocardiograph and other medical electronic devices;

which may be life-threatening. 4. It's not allowed in gas stations, workshops with excessive dusts, gas and other places of flammable and explosive material.

⚠ Warning

- When the high-frequency surgical equipment and the massager are applied on a patient at the same time, the massage pads of the
- massager may cause burns or damage the massager. • When the product is applied close to (1M) a short-wave or microwave treatment device, the output of the massager may be unstable.
- During massage, don't contact wounds and scars with the massage pads. During massage, don't place the massage pads near the heart. Otherwise it will increase the risk of heart fibrillation.
- The product is recommended to be used by one user. In case it's used by two or more users, the massage pads should be cleaned before use. When using the product, the massage pads should be tightly attached
- to the skin in uniform contact. • The massage pads of the massager has the heating function. Hot compress on the same part for a long time (more than 15 minutes) may cause low-temperature burns.
- In case of discomfort during or after usage, stop using it and seek • This product contains a lithium battery. Don't disassemble, modify, immerse the product or throw it into fire, otherwise it may cause
- product damage or explosion. Please consult your doctor for the following situations. Failure to do so may cause an accident or discomfort. 1. Patients who are receiving medical treatment and feel unwell.

2. Patients suffering from abnormal cardiac conditions and brain nerves. 3. Patients whose temperature exceeds 38 degrees (fever).

4. Patients suffering from infectious diseases. 5. Patients subject to skin perceptual disorder or skin abnormalities.

7. Patients with abnormal blood pressure. 8. Patients who cannot express themselves freely.

9. Patients suffering from sensory disorder due to peripheral circulation disorder caused by diabetes, etc.

10. Patients with skin allergies to stainless steel massage pads. Usage is prohibited in the following circumstances

1.Patients with malignant tumors 2. Patients with hemorrhagic tendency 3. Pregnant and lying-in woman

4. Patients with acute suppurative inflammation 5. Pacemaker users

6. Patients with metal implants 7. Near the human heart • Don't apply the product near the heart, above the neck, head, mouth. private parts, or parts with skin diseases and don't massage the soles

of both feet at the same time. In addition, don't compress the internal organs when using. Failure to do so may cause an accident or Don't apply it with other medical devices or use more than two

sets of the product at the same time. Otherwise it may cause nausea or discomfort. Don't apply it for purposes other than massage as described in the user manual. Otherwise accidents and malfunctions may occur.

 Don't use the massage pads with missing or damaged surfaces. Otherwise, accidents or malfunctions may occur. Don't disassemble, repair and modify this product without permission.

⚠ Note

and take out the battery.

In case of discomfort due to product abnormality, stop using it

Otherwise it may cause malfunction or electric shock.

• In case of a strange feeling on your body or skin, consult your doctor and follow your doctor's advice.

08 EN

- During the massage, to move the massage pads or place it on other people, make sure to cut off the power before moving. Otherwise you
- may get a strong electric shock. Don's massage with electronic equipment. Otherwise, it may cause malfunction of electronic devices (such as messy electronic watch time).
- In the event of an abnormal operation, please cut off the power immediately and remove it. Otherwise high temperature or short circuit may occur, leading to an accident. At the same time, contact after-sales service.
- Don't plug or unplug the USB charging cable with wet hands. Otherwise it may cause electric shock or accident.
- Don't allow children or people not capable of expressing consciously to use. Otherwise it may cause accidents or physical discomfort. Don't use the product in a place with high humidity, such as a bathroom,
- or use the product while bathing. Otherwise, you may receive a strong electric shock and damage the product. • Don't use the product when you are in bed. Otherwise, the device may
- malfunction. In addition, the massage pads may be moved to other parts and cause discomfort. Don't use the product when driving. Otherwise it may cause discomfort
- due to excessive stimulation, resulting in an accident. Don't stick the massage pads to the affected area for a long time.
- Otherwise it may cause skin inflammation. When massaging, don't touch the massage pads with the metal parts of the belt, necklace, etc. Otherwise you may get a strong electric shock.
- Don't apply the product on the same part for a long time (more than 30 minutes). Otherwise it may cause muscle fatigue in the massaged • Wear this product correctly. It is used on the cervical spine and is worn
- symmetrically. Otherwise, accidents or discomfort may occur. When disconnecting the USB charging cable from the product, don't pull the cable, but hold the plug to pull. Otherwise it may cause electric
- shock or damage the product. The massage pads of the massager has the heating function. Users who are not sensitive to heat or too sensitive to heat should consult
- a doctor before use. When the product is in use for a long time, the surface will be hot, which is normal.

09 EN

Suggestion

- Make sure that the buttons on the device function normally before use. When the device is used for the first time or not in use for a long time, please charge it according to the user manual to ensure that it can function normally.
- When disposing of this product, please dispose of it in accordance with the relevant environmental protection regulations of the city. Otherwise it may cause environmental pollution.
- This product's target users are people aged 16 or above. Don't use or store this product in the following places:
- Places exposed to direct sunlight Places subject to high temperature and humidity
- Places with water Places with excessive dust Places with smoke or fire

Places subject to vibration Places subject to strong electromagnetism Power control module **EMS** generation

Capacity of battery 1000mAh Rated power 5W

IEC 60335-1 IEC 60335-2-32

K5 Pro

Product and component recycling instructions Recycle instructions Rated voltage 5V === In order to protect the environment, when the product or its parts are discarded, as a consumer you have the responsibility to separate them from the domestic waste and send to a qualified recycling station, by

Name and content of hazardous substances in the product

accordance with relevant national regulations.

which the product or its parts are disassembled and recycled in

art name	Hazardous substances					
	Pb	Hg	Cd	Cr(VI)	PBB	PBDE
lardware	0	0	0	0	0	0
CB board	×	0	0	0	0	0
lastic parts	0	0	0	0	0	0

O: Indicates that the content of the hazardous substance in the homogeneous material of the component is below the limit

requirement specified in GB/T 26572. ×: Indicates that the content of the hazardous substance in at least one homogeneous material of the component exceeds the limit

requirement specified in GB / T 26572. This product has a "green life" of 10 years according to



the product design by the company. The term "green life" is valid only when the product is used under normal Green life | conditions as described in the user manual of this product.

1. All other components and their homogeneous materials, which are not listed in the above table, don't contain harmful substances. 2. All the harmful substances in components and their homogeneous materials listed in the above table are in line with strict requirements of China RoHS on the use of hazardous substances. Please feel assured

1 EN

材质:128g铜版纸 成品尺寸:105×143mm 颜色:Panton 3252C+Coolgray8C+Black 装订方式:骑马钉 编码:208010xxxx

2020-9-26

Note for discarding the product

Please support the recycling of rechargeable battery (lithium ion battery)

The battery used for this product is a dedicated

rechargeable battery. Don't use it on other products. Don't charge the battery that is taken out. Don't throw it into fire and overheat it. Don't pierce it with nail, impact it or disassemble it Don't short-circuit the battery with metal.

 Don't store it with necklace, hairnin. Don't charge, use and place the battery near fire, under the sun and other places of high temperature. The above situation can all cause heating, fire and

Please place the removed rechargeable battery out of the reach of infants and young children to prevent

If swallowed by accident, please go to a doctor

Don't touch the liquid leaking from the battery and

handle it in the following way. When the battery fluid gets into the eyes, it poses danger of blindness. Don't rub your eyes with your hands directly, instead, wash your eyes with water immediately and see a doctor for treatment. When the battery liquid contacts the body or clothing, it may cause skin irritation or injury. Please wash with water and if discomfort persists, you need to go to a doctor

 Disassemble only when the product is discarded. Otherwise, it may cause an accident such as fire, electric shock or injury. If you need to repair (except for replacement of battery), please consult the store where you bought the product.

Lithium-ion batteries used in this product are recyclable and valuable resource. Please hand it over to professional organizations that can carry out rechargeable battery recycling when they are spent and discarded. The battery must be removed from the appliance before it's disposed of; The appliance must be powered off before the removal of battery; The battery shall be disposed of in a safe manner. Note: Use dedicated tools to remove the plastic casing at the non-key side,

official service platform for online customer service, electronic insurance sheet, installation and repair, service schedule, after-sales outlet, membership privilege and latest benefits. Place of origin Warranty period

Please show this card to service outlet or scan the QR code to follow SKG



Purchase date

Inspection date

Signature of repairer

Repair record

Date sent for repair Inspection date Malfunction, repair and handling



Inspector No.: ___